

ESF1 antibody

Product Information

Catalog No.:	FNab02864
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May constitute a novel regulatory system for basal transcription. Negatively regulates ABT1(By similarity).

Immunogen information

Immunogen:	ESF1, nucleolar pre-rRNA processing protein, homolog
Synonyms:	ESF1 homolog ABT1-associated protein ESF1 ABTAP C20orf6
Observed MW:	100-110 kDa
Uniprot ID :	Q9H501

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab02864(ESF1 Antibody) at dilution of 1:2000