

ESRRB antibody

Product Information

Catalog No.:	FNab02871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Nuclear receptor, may regulate ESR1 transcriptional activity. Induces the expression of PERM1 in the skeletal muscle.

Immunogen information


Immunogen:	estrogen-related receptor beta
Synonyms:	Steroid hormone receptor ERR2 ERR beta-2 Estrogen receptor-like 2 Estrogen-related receptor beta (ERR-beta) Nuclear receptor subfamily 3 group B member 2 ESRRB ERRB2 ESRL2 NR3B2
Observed MW:	45-49 kDa
Uniprot ID :	O95718

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-ESRRB (IP:FNab02871, 4ug;
Detection:FNab02871 1:300) with A549 cells
lysate 1040ug.



PC-13 cells were subjected to SDS PAGE
followed by western blot with FNab02871(ESRRB
antibody) at dilution of 1:1000