

## ESRRG antibody

### Product Information

Catalog No.:	FNab02872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Orphan receptor that acts as transcription activator in the absence of bound ligand. Binds specifically to an estrogen response element and activates reporter genes controlled by estrogen response elements(By similarity). Induces the expression of PERM1 in the skeletal muscle.

### Immunogen information

Immunogen:	estrogen-related receptor gamma
Synonyms:	Estrogen-related receptor gamma ERR gamma-2 Estrogen receptor-related protein 3 Nuclear receptor subfamily 3 group B member 3 ESRRG ERR3 ERRG2 KIAA0832 NR3B3
Observed MW:	51 kDa
Uniprot ID :	P62508

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:3000
Image:	



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab02872(ESRRG antibody) at dilution of 1:1000