

ETF1 antibody

Product Information

Catalog No.:	FNab02874
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Directs the termination of nascent peptide synthesis(translation) in response to the termination codons UAA, UAG and UGA. Component of the transient SURF complex which recruits UPF1 to stalled ribosomes in the context of nonsense-mediated decay(NMD) of mRNAs containing premature stop codons.

Immunogen information


Immunogen:	eukaryotic translation termination factor 1
Synonyms:	Eukaryotic peptide chain release factor subunit 1 (Eukaryotic release factor 1, eRF1) Protein C11 TB3-1 ETF1 ERF1 RF1 SUP45L1
Observed MW:	45-48 kDa
Uniprot ID :	P62495

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab02874(ETF1 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02874(ETF1 antibody) at dilution of 1:800