

## ETS2 antibody

### Product Information

Catalog No.:	FNab02879
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor activating transcription. Binds specifically the DNA GGAA/T core motif(Ets-binding site or EBS) in gene promoters and stimulates transcription.

### Immunogen information


Immunogen:	v-ets erythroblastosis virus E26 oncogene homolog 2(avian)
Synonyms:	Protein C-ets-2 ETS2
Observed MW:	53 kDa
Uniprot ID :	P15036

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000
Image:	



IP Result of anti-ETS2 (IP:FNab02879, 3ug;  
Detection:FNab02879 1:600) with HeLa cells  
lysate 5000ug.



K-562 cells were subjected to SDS PAGE  
followed by western blot with FNab02879(ETS2  
antibody) at dilution of 1:1000