

## EVI2A antibody

### Product Information

Catalog No.:	FNab02884
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May complex with itself or/and other proteins within the membrane, to function as part of a cell-surface receptor.

### Immunogen information


Immunogen:	ecotropic viral integration site 2A
Synonyms:	Protein EVI2A Ecotropic viral integration site 2A protein homolog (EVI-2A) EVI2A EVDA EVI2
Observed MW:	26 kDa
Uniprot ID :	P22794

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab02884(EVI2A antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02884(EVI2A antibody) at dilution of 1:300