

EVI2A antibody

Product Information

Catalog No.:	FNab02885
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May complex with itself or/and other proteins within the membrane, to function as part of a cell-surface receptor.

Immunogen information

Immunogen:	ecotropic viral integration site 2a
Synonyms:	Protein EVI2A Ecotropic viral integration site 2A protein homolog (EVI-2A) EVI2A EVDA EVI2
Observed MW:	26 kDa
Uniprot ID :	P22794

Application

Reactivity:	Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02885(Evi2a antibody) at dilution of 1:300