

EXOC7 antibody

Product Information

Catalog No.:	FNab02892
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	exocyst complex component 7
Synonyms:	Exocyst complex component 7 Exocyst complex component Exo70 EXOC7 EXO70 KIAA1067
Observed MW:	74 kDa
Uniprot ID :	Q9UPT5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200; IP: 1:500-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab02892(EXOC7 antibody) at dilution of 1:50



IP Result of anti-Exo70 (IP: FNab02892, 3ug; Detection: FNab02892 1:600) with mouse brain tissue lysate 3000ug.



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02892(EXOC7 antibody) at dilution of 1:500