

## EXPH5 antibody

### Product Information

Catalog No.:	FNab02909
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May act as Rab effector protein and play a role in vesicle trafficking.

### Immunogen information


Immunogen:	exophilin 5
Synonyms:	Exophilin-5 Synaptotagmin-like protein homolog lacking C2 domains b (Slac2-b, Slp homolog lacking C2 domains b) EXPH5 KIAA0624 SLAC2B
Observed MW:	222-240 kDa
Uniprot ID :	Q8NEV8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using FNab02909(EXPH5 Antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab02909(EXPH5 Antibody) at dilution of 1:300