

F2RL3 antibody

Product Information

Catalog No.:	FNab02921
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Receptor for activated thrombin or trypsin coupled to G proteins that stimulate phosphoinositide hydrolysis. May play a role in platelets activation.

Immunogen information


Immunogen:	coagulation factor II(thrombin) receptor-like 3
Synonyms:	Proteinase-activated receptor 4 (PAR-4) Coagulation factor II receptor-like 3 Thrombin receptor-like 3 F2RL3 PAR4
Observed MW:	41 kDa
Uniprot ID :	Q96RI0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse lung tissue slide using FNab02921(F2RL3 Antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab02921(F2RL3 Antibody) at dilution of 1:1000