

F7 antibody

Product Information

Catalog No.:	FNab02924
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Initiates the extrinsic pathway of blood coagulation. Serine protease that circulates in the blood in a zymogen form. Factor VII is converted to factor VIIa by factor Xa, factor XIIa, factor IXa, or thrombin by minor proteolysis. In the presence of tissue factor and calcium ions, factor VIIa then converts factor X to factor Xa by limited proteolysis. Factor VIIa will also convert factor IX to factor IXa in the presence of tissue factor and calcium.

Immunogen information

Immunogen:	coagulation factor VII(serum prothrombin conversion accelerator)
Synonyms:	Coagulation factor VII Proconvertin Serum prothrombin conversion accelerator (SPCA) Factor VII light chain Factor VII heavy chain F7
Observed MW:	49-55 kDa
Uniprot ID :	P08709

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab02924(F7 Antibody) at dilution of 1:300