

## FA2H antibody

### Product Information

Catalog No.:	FNab02926
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Required for alpha-hydroxylation of free fatty acids and the formation of alpha-hydroxylated sphingolipids.

### Immunogen information


Immunogen:	fatty acid 2-hydroxylase
Synonyms:	Fatty acid 2-hydroxylase Fatty acid alpha-hydroxylase Fatty acid hydroxylase domain-containing protein 1 FA2H FAAH FAXDC1
Observed MW:	38-43 kDa
Uniprot ID :	Q7L5A8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab02926(FA2H antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02926(FA2H antibody) at dilution of 1:600