

FABP4 antibody

Product Information

Catalog No.:	FNab02933
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus(By similarity).

Immunogen information


Immunogen:	fatty acid binding protein 4, adipocyte
Synonyms:	Fatty acid-binding protein, adipocyte Adipocyte lipid-binding protein (ALBP) Adipocyte-type fatty acid-binding protein (A-FABP, AFABP) Fatty acid-binding protein 4 FABP4
Observed MW:	15 kDa
Uniprot ID :	P15090

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02933(FABP4 antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab02933(FABP4 antibody) at dilution of 1:300