

FABP7-Specific antibody

Product Information

Catalog No.:	FNab02937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	fatty acid binding protein 7, brain
Synonyms:	Fatty acid-binding protein, brain Brain lipid-binding protein (BLBP) Brain-type fatty acid-binding protein (B-FABP) Fatty acid-binding protein 7 Mammary-derived growth inhibitor related FABP7 BLBP FABPB MRG
Observed MW:	15 kDa
Uniprot ID :	O15540

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02937(FABP7-Specific antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02937(FABP7-Specific antibody) at dilution of 1:500