

## F12 antibody

### Product Information


Catalog No.:	FNab02939
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4E7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	coagulation factor XII(Hageman factor)
Synonyms:	Coagulation factor XII Hageman factor (HAF) Coagulation factor XIIa heavy chain Beta-factor XIIa part 1 Coagulation factor XIIa light chain Alternative names: Beta-factor XIIa part 2 F12
Observed MW:	28 kDa
Uniprot ID :	P00748

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab02939(F12 Antibody) at dilution of 1:200



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02939(F12 antibody) at dilution of 1:1000