

## F13A1 antibody

### Product Information

Catalog No.:	FNab02940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	coagulation factor XIII, A1 polypeptide
Synonyms:	Coagulation factor XIII A chain (Coagulation factor XIIIa) Protein-glutamine gamma-glutamyltransferase A chain Transglutaminase A chain F13A1 F13A
Observed MW:	80-85 kDa
Uniprot ID :	P00488

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human placenta tissue were subjected to SDS  
PAGE followed by western blot with  
FNab02940(Factor XIII A Antibody) at dilution of  
1:4000