

FLAD1 antibody

Product Information

Catalog No.:	FNab02942
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	FAD1 flavin adenine dinucleotide synthetase homolog
Synonyms:	FAD synthase FAD pyrophosphorylase FMN adenylyltransferase Flavin adenine dinucleotide synthase FLAD1
Observed MW:	50 kDa
Uniprot ID :	Q8NFF5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:20-1:200
Image:	

K-562 cells were subjected to SDS PAGE
followed by western blot with FNab02942(FADS
antibody at dilution of 1:600

