

FAF1 antibody

Product Information

Catalog No.:	FNab02946
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9H10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potentiates but cannot initiate FAS-induced apoptosis.

Immunogen information

Immunogen:	Fas(TNFRSF6) associated factor 1
Synonyms:	FAS-associated factor 1 (hFAF1) UBX domain-containing protein 12 UBX domain-containing protein 3A FAF1 UBXD12 UBXN3A
Observed MW:	74kd
Uniprot ID :	Q9UNN5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02946(FAF1 antibody) at dilution of 1:1000