

## FAH antibody

### Product Information

Catalog No.:	FNab02948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


This gene encodes the last enzyme in the tyrosine catabolism pathway. FAH deficiency is associated with Type 1 hereditary tyrosinemia (HT).

### Immunogen information


Immunogen:	fumarylacetoacetate hydrolase (fumarylacetoacetase)
Synonyms:	Fumarylacetoacetase (FAA) Beta-diketonase Fumarylacetoacetate hydrolase FAH
Observed MW:	46 kDa
Uniprot ID :	P16930

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver injury using FNab02948(FAH antibody) at dilution of 1:50



various lysates were subjected to SDS PAGE followed by western blot with FNab02948(FAH antibody) at dilution of 1:1000