

## RFLNA antibody

### Product Information

Catalog No.:	FNab02951
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the regulation of the perinuclear actin network and nuclear shape through interaction with filamins. Plays an essential role in actin cytoskeleton formation in developing cartilaginous cells.

### Immunogen information


Immunogen:	family with sequence similarity 101, member A
Synonyms:	Refilin-A Regulator of filamin protein A (RefilinA) RFLNA FAM101A
Observed MW:	26 kDa
Uniprot ID :	Q6ZTI6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab02951(FAM101A antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02951(FAM101A antibody) at dilution of 1:500