

FAM107A antibody

Product Information

Catalog No.:	FNab02953
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


When transfected into cell lines in which it is not expressed, suppresses cell growth. May play a role in tumor development.

Immunogen information

Immunogen:	family with sequence similarity 107, member A
Synonyms:	Actin-associated protein FAM107A Down-regulated in renal cell carcinoma 1 Protein TU3A FAM107A DRR1 TU3A
Observed MW:	None
Uniprot ID :	O95990

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab02953(FAM107A antibody) at dilution of 1:50