

## FAM111B antibody

### Product Information

Catalog No.:	FNab02954
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	family with sequence similarity 111, member B
Synonyms:	Serine protease FAM111B Cancer-associated nucleoprotein FAM111B CANP
Observed MW:	135 kDa
Uniprot ID :	Q6SJ93

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon slide using FNab02954(FAM111B Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02954(FAM111B antibody) at dilution of 1:500