

RTL8C antibody

Product Information

Catalog No.:	FNab02962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 127, member A
Synonyms:	Retrotransposon Gag-like protein 8C Mammalian retrotransposon derived protein 8C RTL8C CXX1 FAM127A MAR8 MAR8C MART8
Observed MW:	40 kDa
Uniprot ID :	A6ZKI3

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab02962(FAM127A antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02962(FAM127A antibody) at dilution of 1:300