

## NIBAN2 antibody

### Product Information

Catalog No.:	FNab02965
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in apoptosis suppression. May promote melanoma cell invasion in vitro.

### Immunogen information


Immunogen:	family with sequence similarity 129, member B
Synonyms:	Protein Niban 2 Meg-3 Melanoma invasion by ERK (MINERVA) Niban-like protein 1 Protein FAM129B NIBAN2 C9orf88 FAM129B
Observed MW:	100 kDa
Uniprot ID :	Q96TA1

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab02965(FAM129B Antibody) at dilution of 1:50



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab02965(FAM129B Antibody) at dilution of 1:1000