

## EMC9 antibody

### Product Information

Catalog No.:	FNab02969
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	family with sequence similarity 158, member A
Synonyms:	ER membrane protein complex subunit 9 Protein FAM158A EMC9 C14orf122 FAM158A
Observed MW:	23 kDa
Uniprot ID :	Q9Y3B6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab02969(FAM158A Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02969(FAM158A Antibody) at dilution of 1:300