

FHIP2A antibody

Product Information

Catalog No.:	FNab02970
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 160, member B1
Synonyms:	FHF complex subunit HOOK interacting protein 2A (FHIP2A) FHIP2A FAM160B1 KIAA1600
Observed MW:	105 kDa
Uniprot ID :	Q5W0V3

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-FAM160B1 (IP:FNab02970, 5ug;
Detection:FNab02970 1:300) with PC-3 cells
lysate 1500ug.



PC-3 cells were subjected to SDS PAGE followed
by western blot with FNab02970(FAM160B1
Antibody) at dilution of 1:600