

FAM3D antibody

Product Information

Catalog No.:	FNab02983
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	family with sequence similarity 3, member D
Synonyms:	Protein FAM3D FAM3D
Observed MW:	30 kDa
Uniprot ID :	Q96BQ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab02983(FAM3D antibody) at dilution of 1:6000