

## STRIP2 antibody

### Product Information

Catalog No.:	FNab02984
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	family with sequence similarity 40, member B
Synonyms:	Striatin-interacting protein 2 Protein FAM40B STRIP2 FAM40B KIAA1170
Observed MW:	95 kDa
Uniprot ID :	Q9ULQ0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab02984( FAM40B Antibody) at dilution of 1:600