

CYRIB antibody

Product Information

Catalog No.:	FNab02987
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 49, member B
Synonyms:	CYFIP-related Rac1 interactor B L1 CYRIB CYRI FAM49B
Observed MW:	37 kDa
Uniprot ID :	Q9NUQ9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab02987(FAM49B antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02987(FAM49B antibody) at dilution of 1:500