

GARIN5B antibody

Product Information

Catalog No.:	FNab02992
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 71, member E2
Synonyms:	Golgi-associated RAB2 interactor protein 5B GARIN5B C19orf16 FAM71E2
Observed MW:	34 kDa
Uniprot ID :	Q8N5Q1

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab02992(FAM71E2 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02992(FAM71E2 antibody) at dilution of 1:1000