

## GARIN1A antibody

### Product Information

Catalog No.:	FNab02993
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	family with sequence similarity 71, member F2
Synonyms:	Golgi-associated RAB2 interactor protein 1A GARIN1A FAM137B FAM71F2
Observed MW:	40-45 kDa
Uniprot ID :	Q6NXP2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab02993( FAM71F2 Antibody) at dilution of 1:2000