

LRATD2 antibody

Product Information

Catalog No.:	FNab02997
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 84, member B
Synonyms:	Protein LRATD2 Breast cancer membrane protein 101 LRAT domain-containing 2 Protein FAM84B Protein NSE2 LRATD2 BCMP101 FAM84B NSE2
Observed MW:	40 kDa
Uniprot ID :	Q96KN1

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab02997(FAM84B Antibody) at dilution of 1:200



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with FNab02997(FAM84B Antibody) at dilution of 1:600