

FAM8A1 antibody

Product Information

Catalog No.:	FNab02998
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	family with sequence similarity 8, member A1
Synonyms:	Protein FAM8A1 Autosomal highly conserved protein FAM8A1 AHCP
Observed MW:	44 kDa
Uniprot ID :	Q9UBU6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab02998(FAM8A1 Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02998(FAM8A1 Antibody) at dilution of 1:1000