

FANCA-Specific antibody

Product Information

Catalog No.:	FNab03004
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Fanconi anemia, complementation group A
Synonyms:	Fanconi anemia group A protein (Protein FACA) FANCA FAA FACA FANCH
Observed MW:	163 kDa
Uniprot ID :	O15360

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03004(FANCA-Specific antibody) at dilution of 1:500