

FANCG antibody

Product Information

Catalog No.:	FNab03008
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. May be implicated in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability. Candidate tumor suppressor gene.

Immunogen information

Immunogen:	Fanconi anemia, complementation group G
Synonyms:	Fanconi anemia group G protein (Protein FACG) DNA repair protein XRCC9 FANCG XRCC9
Observed MW:	65-68 kDa
Uniprot ID :	O15287

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03008(FANCG antibody) at dilution of 1:100