

FARSA antibody

Product Information

Catalog No.:	FNab03014
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	phenylalanyl-tRNA synthetase, alpha subunit
Synonyms:	Phenylalanine--tRNA ligase alpha subunit CML33 Phenylalanyl-tRNA synthetase alpha subunit (PheRS) FARSA FARS FARSL FARSLA
Observed MW:	57 kDa
Uniprot ID :	Q9Y285

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab03014(FARSA antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03014(FARSA antibody) at dilution of 1:500