

FSCN1 antibody

Product Information

Catalog No.:	FNab03018
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	fascin homolog 1, actin-bundling protein(Strongylocentrotus purpuratus)
Synonyms:	Fascin 55 kDa actin-bundling protein Singed-like protein p55 FSCN1 FAN1 HSN SNL
Observed MW:	50-55 kDa
Uniprot ID :	Q16658

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab03018(Fascin Antibody) at dilution of 1:200



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab03018(Fascin Antibody) at dilution of 1:600