

FASTKD2 antibody

Product Information

Catalog No.:	FNab03020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays an important role in assembly of the mitochondrial large ribosomal subunit(PubMed:25683715).

Immunogen information


Immunogen:	FAST kinase domains 2
Synonyms:	FAST kinase domain-containing protein 2, mitochondrial FASTKD2 KIAA0971
Observed MW:	65-75 kDa
Uniprot ID :	Q9NYY8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-FASTKD2 (IP:FNab03020, 4ug;
Detection:FNab03020 1:300) with mouse heart
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab03020(FASTKD2 antibody) at dilution of
1:1000