

FATE1 antibody

Product Information

Catalog No.:	FNab03022
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	fetal and adult testis expressed 1
Synonyms:	Fetal and adult testis-expressed transcript protein Cancer/testis antigen 43 (CT43) Tumor antigen BJ-HCC-2 FATE1 FATE
Observed MW:	25 kDa
Uniprot ID :	Q969F0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab03022(FATE1 Antibody) at dilution of 1:800