

FAU antibody

Product Information

Catalog No.:	FNab03024
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Finkel-Biskis-Reilly murine sarcoma virus(FBR-MuSV) ubiquitously expressed
Synonyms:	Ubiquitin-like FUBI-ribosomal protein eS30 fusion protein FAU ubiquitin like and ribosomal protein S30 fusion Ubiquitin-like protein FUBI Small ribosomal subunit protein eS30 Alternative names: 40S ribosomal protein S30 FAU
Observed MW:	14 kDa
Uniprot ID :	P35544

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab03024(FAU Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03024(FAU antibody) at dilution of 1:300