

## FBF1 antibody

### Product Information

Catalog No.:	FNab03026
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Keratin-binding protein required for epithelial cell polarization. Involved in apical junction complex(AJC) assembly via its interaction with PARD3. Required for ciliogenesis.

### Immunogen information


Immunogen:	Fas(TNFRSF6) binding factor 1
Synonyms:	Fas-binding factor 1 (FBF-1) Protein albatross FBF1 ALB KIAA1863
Observed MW:	69-70 kDa
Uniprot ID :	Q8TES7

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03026(FBF1 antibody) at dilution of 1:50



IP Result of anti-FBF1 (IP:FNab03026, 4ug; Detection:FNab03026 1:300) with mouse liver tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03026(FBF1 antibody) at dilution of 1:400