

## FBXL13 antibody

### Product Information

Catalog No.:	FNab03032
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Substrate-recognition component of the SCF(SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

### Immunogen information


Immunogen:	F-box and leucine-rich repeat protein 13
Synonyms:	F-box and leucine-rich repeat protein 13 Dynein regulatory complex subunit 6 F-box/LRR-repeat protein 13 FBXL13 DRC6 FBL13
Observed MW:	120 kDa
Uniprot ID :	Q8NEE6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab03032(FBXL13 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03032(FBXL13 antibody) at dilution of 1:500