

FBXL18 antibody

Product Information

Catalog No.:	FNab03036
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Substrate-recognition component of the SCF(SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Immunogen information

Immunogen:	F-box and leucine-rich repeat protein 18
Synonyms:	F-box/LRR-repeat protein 18 F-box and leucine-rich repeat protein 18 FBXL18 FBL18
Observed MW:	None
Uniprot ID :	Q96ME1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab03036(FBXL18 Antibody) at dilution of 1:50