

## FBXO3 antibody

### Product Information

Catalog No.:	FNab03043
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Substrate recognition component of the SCF(SKPI-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Mediates the ubiquitination of HIPK2 and probably that of EP300, leading to rapid degradation by the proteasome. In the presence of PML, HIPK2 ubiquitination still occurs, but degradation is prevented. PML, HIPK2 and FBXO3 may act synergically to activate p53/TP53-dependent transactivation.

### Immunogen information

Immunogen:	F-box protein 3
Synonyms:	F-box only protein 3 FBXO3 FBX3
Observed MW:	48-50 kDa
Uniprot ID :	Q9UK99

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03043(FBOXO3 antibody) at dilution of 1:800