

FBXO5 antibody

Product Information

Catalog No.:	FNab03047
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Regulates progression through early mitosis by inhibiting the anaphase promoting complex/cyclosome(APC). Binds to the APC activators CDC20 and FZR1/CDH1 to prevent APC activation. Can also bind directly to the APC to inhibit substrate-binding.

Immunogen information


Immunogen:	F-box protein 5
Synonyms:	F-box only protein 5 Early mitotic inhibitor 1 FBXO5 EMI1 FBX5
Observed MW:	56 kDa, 45 kDa
Uniprot ID :	Q9UKT4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab03047(FBXO5 Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03047(FBXO5 antibody) at dilution of 1:800