

## APBB1 antibody

### Product Information

Catalog No.:	FNab03069
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	amyloid beta(A4) precursor protein-binding, family B, member 1(Fe65)
Synonyms:	Amyloid beta precursor protein binding family B member 1 Amyloid-beta A4 precursor protein-binding family B member 1 Protein Fe65 APBB1 FE65 RIR
Observed MW:	77 kDa
Uniprot ID :	O00213

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03069(FE65 antibody) at dilution of 1:300