

FEM1B antibody

Product Information

Catalog No.:	FNab03072
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit. Involved in apoptosis by acting as a death receptor-associated protein that mediates apoptosis. Also involved in glucose homeostasis in pancreatic islet. Functions as an adapter/mediator in replication stress-induced signaling that leads to the activation of CHEK1.

Immunogen information


Immunogen:	fem-1 homolog b
Synonyms:	Protein fem-1 homolog B (FEM1b) FEM1-beta Fem-1-like death receptor-binding protein alpha Fem-1-like in apoptotic pathway protein alpha (F1A-alpha) FEM1B F1AA KIAA0396
Observed MW:	70 kDa
Uniprot ID :	Q9UK73

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab03072(FEM1B antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03072(FEM1B antibody) at dilution of 1:1000