

## FEM1C antibody

### Product Information

Catalog No.:	FNab03074
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit.

### Immunogen information


Immunogen:	fem-1 homolog c(C. elegans)
Synonyms:	Protein fem-1 homolog C (FEM1c) FEM1-gamma FEM1C KIAA1785
Observed MW:	69 kDa
Uniprot ID :	Q96JP0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab03074(FEM1C Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03074(FEM1C antibody) at dilution of 1:1000