

AHSG antibody

Product Information

Catalog No.:	FNab03080
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	alpha-2-HS-glycoprotein
Synonyms:	Alpha-2-HS-glycoprotein Alpha-2-Z-globulin Ba-alpha-2-glycoprotein Fetuin-A Alpha-2-HS-glycoprotein chain A Alpha-2-HS-glycoprotein chain B AHSG FETUA
Observed MW:	60 kDa
Uniprot ID :	P02765

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab03080(AHSG Antibody) at dilution of 1:200



human liver tissue were subjected to SDS PAGE followed by western blot with FNab03080(AHSG antibody) at dilution of 1:500