

FEV antibody

Product Information

Catalog No.:	FNab03083
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a transcriptional regulator. According to PubMed:12761502, it functions as a transcriptional repressor. Functions in the differentiation and the maintenance of the central serotonergic neurons. May play a role in cell growth.

Immunogen information

Immunogen:	FEV(ETS oncogene family)
Synonyms:	Protein FEV Fifth Ewing variant protein PC12 ETS domain-containing transcription factor 1 (PC12 ETS factor 1, Pet-1) FEV PET1
Observed MW:	25 kDa
Uniprot ID :	Q99581

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



fetal human brain tissue were subjected to SDS
PAGE followed by western blot with
FNab03083(FEV Antibody) at dilution of 1:1000